

Foreclip Device Instructions

- The Foreclip **can be worn for ~1 hour**. The device is recommended for home usage so you can readjust. The device is for skin below the scar line only
- **Grip the clips onto the skin without them being connected**. Once they are on then connect them and set tension. If there is pain after the clips are connected - it means the tension is too high
- when putting on the device, make sure you **grip a significant amount of skin (skin should go to top of cushions)** this is how you can generate more tension. Do not clip directly onto any veins.
- **Top screws can be adjusted to change grip strength**
- To change grippers use a flathead screwdriver to pry them out. Eventually, grip cushions will need to be replaced

Common Issues

- **Tension set too high**. If the tension is set too high the device will be painful. To properly discern the pinching feeling (where device grips skin) from the stretching feeling (tensioned skin), clamp the device onto your skin but do not connect the clips.
- **Not enough skin under the gripper**. The grippers are simulating what your fingers would typically do for Manual Method 1. You must grip a significant amount of skin to generate tension. If this is incorrect, the device is ineffective.

Disclaimer

The safety of this product is not guaranteed. Users must be careful to not injure themselves with the device, if you feel any pain - stop immediately. Do not set the tension too high as that can result in injury. NOT FOR USE ON INNER SKIN (i.e. skin above scar line). The clips may feel a bit strange at first but after a minute it becomes more comfortable